

This Pandemic Season, Trust RONCO for your PPE needs

BLURITE XPL

Nitrile Examination Gloves

- · Protection against bio-hazardous and infectious material
- 100% powder-free to eliminate risks of cross contamination
- 12" long for extended forearm protection







Packaging:

50 gloves x 20 dispensers = 1000 dispensers per case

SOL-FIT

Nitrile Reusable Glove, Unlined

- Superior chemical resistance, ideal for cleaning and disinfecting operations
- Excellent abrasion resistance
- Wet / Dry / Light oil grip



Item #: 29-903-XX Sizes: Small to X-Large

Packaging:

12 pairs x 12 bags = 144 pairs per case

TYCHEM®

QC Coveralls*

- Storm flap covers zipper that can be sealed with adhesive strip to prevent intrusion
- . Compliant for North American standards for blood-borne pathogens (ASTM F1670 and ASTM F1671)



Item #: DP-QC122BYL-XX Sizes: Medium to 5X-Large

Packaging: 4 coveralls per case

FLAT FOLD SURGICAL N95 RESPIRATOR

- · Individually packaged to keep them from exposure, soiling and contamination
- Offers > 99.9% Bacterial Filtration Efficiency (BFE) and > 99.5% Particle Filtration Efficiency (PFE)
- . Conforms to ASTM F1862 requirements for fluid penetration and splash resistance



Item #: SP-14110451

Packaging: Individually wrapped, 200 per case

A600 Splash Goggles*

12" long

mil thick

Examination Grade


Quality

- · Protection against chemical splash, airborne particulates and bio-hazardous material
- · Anti-fog coating
- . Meets ANSI Z87 standards and certified to the requirements of CSA Z94.3 standards



Item #: SP-A610S One size, Indirect ventilation

Packaging: 20 goggles per case



This is **RONCO's** suggested product offering for infection control. For specific enquiries on the right product for the application, please contact **RONCO.**

> Tel: 905.660.6700 Email: ronco@ronco.ca www.ronco.ca



